

We do Beds



	ALEGIO NG	CS8 FLOORLINE SINGLE	CS8 FLOORLINE KING SINGLE	CS7	CS7 KING SINGLE
Bed Length:	2110mm	2110mm	2110mm	2140mm	2140mm
Bed Width	900mm	920mm	1090mm	900mm	1090mm
Adjustable Height:	280 - 800mm	150 - 680mm	150 - 680mm	170 - 760mm	170 - 760mm
Max. User Weight:	135kg	250kg	250kg	226kg	226kg
Mattress Length:	1970mm	1970mm	1970mm	1970mm	1970mm
Mattress Width:	880mm	880mm	1070mm	880mm	1070mm
Bed Weight	78kg	119kg	125kg	86kg	92kg



	ETUDE PLUS	ETUDE PLUS LOW	OCTAVE	1600	SCANBETA NG
Bed Length:	2190mm	2190mm	2221mm	2060mm	1680 or 1880mm
Bed Width	1020mm	1020mm	1190mm	920mm	780 or 880mm
Adjustable Height:	330 - 730mm or 400 - 800mm	200 - 600mm or 270 - 670mm	400 - 800mm	200 - 670mm	300 - 900mm
Max. User Weight:	165kg	165kg	350kg	165kg	70kg
Mattress Length:	1970mm	1970mm	2000mm	1970mm	1600 or 1800mm
Mattress Width:	880mm	880mm	1160mm	880mm	700 or 800mm
Bed Weight	77kg	77kg	170kg	170kg	91kg or 99.5kg

For more information on the Invacare Beds range call **1800 460 460** or visit www.invacare.com.au



Yes, you can.®

Making Life's Experiences Possible