



THE pegasus PRO

COMFORT, SAFETY & RELIABILITY
WITH ADDED PERFORMANCE

With enhanced battery performance, superior suspension and large 12" wheels, the Pegasus PRO is ideal for users who need that extra performance.



Yes, you can.®




THE pegasus PRO

COMFORT, SAFETY & RELIABILITY
WITH ADDED PERFORMANCE

Features and benefits

- The advanced suspension system brings better performance and gives a smoother ride over a variety of surfaces and terrains.
- Features a powerful motor and 12" pneumatic tyres to ensure you get where you need to go in comfort and style.
- An on-board electronic speed reduction system automatically reduces speed around bends.
- A two-step disengaging lever prevents the scooter from free-wheeling if it is knocked out of drive mode.
- The electronics and motor have additional protection to guard against water and corrosion, safeguarding the vital workings of the scooter.

Tech Specs



Driving unit width
650 mm



Total length
1320 mm



Max driving distance
52 km



Max driving speed
10 km/h

Maximum driving distance is dependent on driving habits, terrain and battery age.

For further information and full technical specifications please see your dealer or visit www.invacare.co.nz or www.invacare.com.au



Yes, you can.™

