



THE pegasus METRO

THE PERFECT MIX OF COMFORT,
SAFETY AND RELIABILITY

The ergonomically designed Pegasus METRO has been designed from the ground up to focus on key aspects important to both new and seasoned scooter users.



Yes, you can.®




THE pegasus METRO

THE PERFECT MIX OF COMFORT,
SAFETY AND RELIABILITY

Features and benefits

- The suspension system has been designed to navigate a variety of everyday surfaces and terrains.
- In addition to a unique lighting system, the Pegasus METRO is equipped with an electronic speed reduction feature that automatically reduces speed around bends.
- A two-step disengaging lever prevents the scooter from free-wheeling if it is knocked out of drive mode.
- Electronics and the motor are protected from corrosion and water which safeguard the vital workings of the scooter against dirt preventing any unnecessary malfunctions.

Tech Specs



Driving unit width
660 mm



Total length
1270 mm



Max driving distance
54 km



Max driving speed
10 km/h

Maximum driving distance is dependent on driving habits, terrain and battery age.

For further information and full technical specifications please see your dealer or visit www.invacare.co.nz or www.invacare.com.au

