

# Action 2NG

SIMPLE, ROBUST AND EFFICIENT



The Action 2NG provides easy handling making it the wheelchair ideal for basic needs with quality components to give effortless rolling performance both in doors and out. A great bonus to the Action 2NG is it's cross compatibility throughout the Action range, offering familiarity and flexibility to a hugely successful Action family.



Yes, you can.®

Making Life's Experiences Possible



# Action 2NG

With a total weight starting from 15kg, the Invacare Action 2NG fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle, and being lightweight makes transportation trouble-free.

Cross-compatibility between a number of components on the Action range make this an ideal fleet product that can be easily adapted to meet individual needs. Available in self-propelled and transit versions.

## Features



Decrease seat depth by 50mm for customised seat size.



The 80° footrest swings in and out for easy transfers in or out of the chair.



The Axle position can be adjusted in height by 5cm.



Lightweight foldable frame.

## Tech Specs



**Seat width**  
40.5cm, 45.5cm,  
48cm



**Seat depth**  
40-43cm



**Seat height**  
48.5cm



**Backrest height**  
43cm



**Total width**  
seat width + 20cm  
(self-propel)  
seat width + 14.5cm  
(transit)



**Total length**  
108cm



**Max. user weight**  
125kg



**Transport weight**  
12.1kg (transit)  
8.9kg (self-propel)



**Total weight unloaded**  
14.9kg (transit)  
16.2kg (self-propel)



**Turning radius**  
87cm



**Max. safe slope**  
9°



**Crash Test Approved**  
Fulfills the requirements of ISO 7176-19

FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email sales@invacare.co.nz or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email orders@invacare.com.au or visit [www.invacare.com.au](http://www.invacare.com.au)

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2024



Yes, you can.®