

Rea Azalea

EXCELLENT SEATING COMFORT



The Rea Azalea is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the seat position to be altered, ensuring maximum comfort at all times. This combined with good driving characteristics makes Azalea the ideal choice.



Yes, you can.®



Confidence in posture solutions



Rea Azalea

The Rea Azalea has seat tilt and a 30° backrest recline. Both tilt and recline can be easily adjusted.

Invacare's Dual Stability System (DSS) helps to prevent Decubitus thanks to its non-friction body positioning, which includes a biomechanical positioning backrest hinge that follows the user's natural back movement when adjusting the backrest angle.

Features



Flex 3, backrest plate. Apart from being height adjustable it is now width adjustable as well.



Option for electric adjustment of tilt and recline via a hand control that can be used by the attendant or user.

Tech Specs



Seat width
39 - 54cm



Seat depth
43 - 50cm



Seat height
45cm



Total width
seat width + 23cm



Total length
102cm



Legrest angle
Adjustable



Tilt Range
25°



Max. user weight
135kg



Transport weight
16.5kg



Total weight unloaded
34kg



Turning Radius
77cm



Max. safe slope
7°



Crash Test Approved
Fulfills the requirements of ISO 7176-19

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 07/2024



Yes, you can.®