


EBL Xtra Portable Ramps

DEPLOYED IN SECONDS



The **Invacare EBL Xtra portable ramps** are a new range that have a safe working load of 400 kg and a handy 2 fold design. Based on the EBL range the **Invacare EBL Xtra portable ramps** have a 5 cm Edge Barrier design for safety and security to prevent steering off the edge of the ramp.



Yes, you can.®

Making Life's
Experiences Possible


EBL Xtra Portable Ramps



- **Portable** - Invacare Portable Ramps can be quickly folded into a single compact unit - allowing for easy storage or transport in any vehicle or home environment
- **Lightweight** - Made from a unique combination of uni-directional glass and carbon fibre materials. Invacare Portable Ramps are very strong, yet light enough to be portable
- **Compact** - Dual panel design held together with a patented hinge system connected by webbing enabling Invacare Portable Ramps to fold with a bi-fold action, creating a smaller width
- **Safe** - Clearly marked and grip tape covered, Invacare Portable Ramps are also fitted with pads at the top and base. Invacare Portable Ramps provide a non-slip surface, for all weather conditions
- **Strong** - Invacare Portable Ramps are tested to AS/NZ Standards, up to 400kg!* for security and peace of mind
- **Simple to use** - Invacare Portable Ramps require no assembly allowing the ramp to be deployed in seconds, and immediately ready for use!

Features

- Marked 'High Side' and 'Low Side' to show head and foot of ramp
- Hinge system the ramp to be folded to a single compact unit
- Ramp material is non-conductive tested according to AS1418.10 - 2011.
- Grip tape and pads at both ends
- Four foot-straps positioned at corners for operator to stabilize the ramp in use



Tech Specs

	EBL Xtra 120 DX120-W05L43N51-M000	EBL Xtra 150 DX150-W05L43N52-M000
Length	120cm	150cm
Width	89.5cm	89.5cm
Height	6.8cm	6.8cm
Weight	6.7kg	8.3kg
Safe working load tested according to AS/NZS 3856.1:1998(R2016)	400kg	400kg

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz
Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 04/2024



Yes, you can.®