

# Mini Back

DESIGNED WITH KIDS IN MIND



The Invacare® Matrix® Mini Back features Mini Set hardware, a miniature version of the popular Easy Set hardware, designed with kids in mind. The Mini back is contoured to fit kids of all shapes and sizes - combining exceptional support in an appealing low profile package.



**matrix**®  
SEATING SERIES

# Mini Back



## Cover:

- Foam and cover overlap the shell for comfort and protection.
- Water resistant and breathable, with additional options of bright pink and blue
- Includes removable lumbar positioning pad

## Back shell:

- Light weight, durable aluminum shell, designed for easy attachment of optional lateral supports and headrest.

## Hardware:

- “Mini Set” hardware installs and adjusts in minutes. Small attachment point for minimal interference on back canes.

## Accessories

Fixed Lateral pads



Mini Swing-away Laterals



Mini-Set Mounting Hardware



Chest Harness Interface



Cover Colours Available










Pink Lollipop Startex



Blizzard Blue Startex

## Tech Specs

<b>Mini Back</b>	 <p><b>Height</b> 7", 10", 12", 14"*</p>	 <p><b>Width</b> 10" - 14"</p>	 <p><b>Depth adjustment</b> 42" forwards, 2" backwards</p>	 <p><b>Angle adjustment</b> 40° range</p>	 <p><b>Product weight</b> from 600g (1.3lb)</p>	 <p><b>Weight Capacity</b> 68kg (150lb)</p>	 <p><b>Limited Warranty</b> 2 years</p>
------------------	---	---	---	--	--	--	--

### FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [orders@invacare.com.au](mailto:orders@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2024



**Yes, you can.®**