

ROVI A3

EXCLUSIVELY WITH

MPS
maxx



STAND UP FOR YOURSELF

The all **NEW MPS Maxx - Modular Power Standing System** is available exclusively on the **ROVI A3** base. Our innovative, highly adjustable **standUP Maxx** system combines a power standing function with a full range of power positioning options, all on a **mid-wheel drive** base that delivers a uniquely compact footprint. An outstanding combination of **independence, function and accessibility.**



Yes, you can.®



ADVANCING MOBILITY.



Accessibility within reach

In addition to the enhanced mobility and freedom offered by our MPS Maxx system, the medical benefits of adjusting to a standing position throughout the day are numerous, including better bone health, greater joint mobility, pain management, enhanced renal/bowel function and decreased risk for skin tissue breakdown.



7" Elevating Seat



45° CG Tilt



165° Recline



Power Centre Mount

ROVI A3 Base Specifications

Suspension System	Active Ride Control (ARC)
Base Width	25.5" (65cm)
Base Length	36.5" (93cm)
Drive Wheel Size	14" (35.5cm)
Caster Wheel Size (front/rear)	6" (15cm)
Base only Turning Radius	20.5" (52cm)
Drivetrain	4-pole EAD Motors
Batteries	2 x M34
Battery Charger	8A, Stand Alone
Range ⁽¹⁾	28km (17.7 miles)
Base Weight (inc batteries)	107kg (236lb)
Ground Clearance	3" (7.5cm)
Maximum Speed w/MPS System	9.3kmh (5.8mph)

⁽¹⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

MPS System Features

CG Tilt	45°
Recline with ESR*	165°
Standing	✓
Power Adjustable Seat Height	7" (18cm)
Power Elevating Belt Drive Center Mount Articulated Foot Platform	✓
Seat Width	16"-20" (40.5-50.5cm)
Seat Depth	16"-20" (40.5-50.5cm)
Seat-to-Floor	As low as 18.5"
Weight Capacity	113kg (250lb)
MPS System Weight ⁽¹⁾	95kg (210lb)

*ESR Extended Shear Reduction

(1) Actual weight will vary based on seating options selected.

Colour Options



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 04/2024



Yes, you can.®