



**TRANSACTION TYPE:**

Order       Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Post Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

Code	Description	Price
<input type="radio"/> BQ1529487-AF	Action2 NG, 40.5cm Seat Width, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/> BQ1529466-AF	Action2 NG, 40.5cm Seat Width, Self-Propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<input type="radio"/> BQ1529491-AF	Action2 NG, 45.5cm Seat Width, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/> BQ1529470-AF	Action2 NG, 45.5cm Seat Width, Self-Propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<input type="radio"/> BQ1529492-AF	Action2 NG, 48cm Seat Width, Transit <i>Equipped with 12" transit wheels inc. fixed axles</i>	
<input type="radio"/> BQ1529485-AF	Action2 NG, 48cm Seat Width, Self-Propelled <i>Equipped with 24" self propelling wheels inc. quick-release axles</i>	
<b>FRAME INCLUSIONS</b>		
CZA 0980	Anthracite Black colour	N/C
CZA 0996	Outside Alu fork, fixed angle caster mount	N/C
CZA 1010	200x32 (8" x 1 1/4") castors, solid tyres	N/C
CZA	Solid profiled rear tyres	N/C
CZA 1270	Push to lock brakes	N/C
<b>SEAT INCLUSIONS</b>		
CZA 0253	Depth adjustable (40 to 43cm) seat upholstery with velcro and handle	N/C
CZA 0270	48.5cm seat height	N/C
<b>BACK INCLUSIONS</b>		
CZA 0455	Folding backrest <i>2 heights 40/43cm (delivered at 43cm)</i>	N/C
CZA 0470	Standard padded nylon sling with flap	N/C
<b>LEGREST INCLUSIONS</b>		
CZA 0510	80° straight swing in/out hanger RH	N/C
CZA 0560	80° straight swing in/out hanger LH	N/C
CZA 0600	Flip-up footplate RH	N/C
CZA 0620	Flip-up footplate LH	N/C
<b>ARMREST INCLUSIONS</b>		
CZA 0650	Removable full-length armrest RH <i>Adjustable 20/23/26cm height increments, inc. locking pin</i>	N/C
CZA 0750	Removable full-length armrest LH <i>Adjustable 20/23/26cm height increments, inc. locking pin</i>	N/C
<b>SAFETY INCLUSIONS</b>		
CZA 1355	Buckle belt, mounted	N/C
CZA 1370	Rear anti-tippers, pair	N/C
CZA 1385	Rear tie-down hooks	N/C
<b>OPTIONAL EXTRAS</b>		
<b>ELEVATING LEGRESTS</b>		
<input type="radio"/>	Angle adj legrest with angle/depth adj footplate RH <input type="radio"/> 40cm (AP1543228R) <input type="radio"/> 45cm (AP1543230R) <input type="radio"/> 48cm (AP1543231R)	
<input type="radio"/>	Angle adj legrest with angle/depth adj footplate LH <input type="radio"/> 40cm (AP1543228L) <input type="radio"/> 45cm (AP1543230L) <input type="radio"/> 48cm (AP1543231L)	
<input type="radio"/>	Amputee legrest RH (1442990)	
<input type="radio"/>	Amputee legrest LH (1442991)	
<b>TRAY</b>		
<input type="radio"/>	Tray table <input type="radio"/> 40cm (1450922) <input type="radio"/> 45cm (1450924) <input type="radio"/> 48cm (1450925)	
<b>HEADREST</b>		
<input type="radio"/>	Neckrest inc. backbrace <input type="radio"/> 40cm (AP1554576) <input type="radio"/> 45cm (AP1554578) <input type="radio"/> 48cm (AP1554579)	
<b>CANE HOLDER</b>		
<input type="radio"/>	Cane holder (0572-02056)	



Print Form

Email Form

Clear Form

ISO, MDA, TUV Logos  
 The product includes all standard features  
 Maximum user weight: 125kg  
 Crossbar Folding Mechanism, Seat Angle Possibility: Delivered Always 3°  
 Rear Wheel Position : Standard, Passive (Delivered Always In Standard)  
 Revised July 2014

