

# KSL

A LITTLE OUT OF THE ORDINARY



The KSL is the full küschall experience with the lightest and most dynamic wheelchair we have ever built. To do this, we had to think about every element. Not just to make it lighter and sturdier - but better as well. The result is something more than a wheelchair. It is The KSL.



Yes, you can.®

*küschall*®



# KSL

In creating this wheelchair, we have fused state-of-the-art technology and design to create a highly dynamic rigid wheelchair for independent active users. With the name The KSL, Kuschall combines the K, a symbol for kuschall® quality and design with “SL”, for Super Light.

The Kuschall-Team have made the sleek and dynamic design of The KSL priority, not only keeping unseemly parts to a minimal but reducing the total weight of the chair to in its standard configuration to less than 7kg.

## Features



Customizable bearing block cap with text or graphic.



Minimal frame design.



Beautiful details with swiss precision.



Carbon camber bar.

## Tech Specs



**Seat width**  
34 - 44cm



**Seat depth**  
38.5 - 48.5cm



**Seat height**  
front: 45 - 53cm  
rear: 40 - 48cm



**Backrest height**  
27 - 42cm



**Frame angle**  
75° / 90°  
V-frame or tapered



**Total width**  
seat width + 17cm



**Total length**  
75°: 80cm  
90°: 73cm



**Max. user weight**  
100 kg



**Transport weight**  
w/out rear wheels  
4.3 kg



**Total weight unloaded**  
6.8 kg



**Turning radius**  
92cm



**Max. safe slope**  
7°



FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [sales@invacare.com.au](mailto:sales@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

**Asia:** Phone +61 2 8839 5330, email [exportsales@invacare.com.au](mailto:exportsales@invacare.com.au) or visit [www.invacare.com](http://www.invacare.com)

**Thailand:** Phone +66 0 2108 8555 or visit [www.invacare.com](http://www.invacare.com)



Yes, you can.®